

PROJECT DEVELOPMENT PHASE

SPRINT 2

Development of Web User Interface in Node-RED platform

Date	05 November 2022
Team ID	PNT2022TMID18885
Project Name	Project - Personal Assistance for Seniors Who Are Self-Reliant

Sprint	Functional Requirement (Epic)	User Story Number	User Story / Task	Story Points	Priority	Team Members
Sprint-1	Registration: Creation of IBM services like NodeRED, CloudantDB, TTS Service and design of IoT system	USN-1	As a user, I should login into my IBM Cloud account.	2	High	Anne Angelina J, Kawin M
Sprint-2	Web UI: Creating web UI using node-red and connect it to IBM Cloudant db	USN-2	As a user, I should be able to feed the medicine name and intake time in the web UI	2	High	Akshayasri S, Bhavani R K
Sprint-3	Hardware implementation: Developing Python code to retrieve data from cloudant db to send that data to IoT device at the appropriate time	USN-3	As a user, I should be able to send the medicine name to the IoT device at the scheduled time	2	High	Akshayasri S, Kawin M
Sprint-4	Software implementation: Converting the data received from cloud as voice using IBM Text to Speech service	USN-4	As a user, I must be able to hear the medicine name which is to be taken at the appropriate time	2	High	Anne Angelina J, Bhavani R K

Objectives:

- i) To create a form UI in Node-Red platform to enter the medicine name and time of intake.
- ii) To send the medicine name at the scheduled time.
- iii) Complete Web UI flow
- iv) Deploying Web UI

i) Creating a form in Node-Red to enter medicine details:

- Web UI can be created by using form UI from Node-Red Dashboard
- The Dashboard can be installed from node-red palette
- The medicine name, intake time and date is entered in the form

- The medicine details are wired and stored to the IBM cloudant db database with the help of a function.

Name

Name

SetupOn StartOn MessageOn Stop

```
1 var d=msg.payload.date
2 var t=msg.payload.time
3 msg.payload= {
4   "_id":d+" "+t,
5   "medicine":msg.payload.medicine,
6 }
7 return msg;
8
9
```

- The below function is used to reset or clear the form after a delay to add more medicine details

Edit function node

DeleteCancelDone

Properties


Name

Name


SetupOn StartOn MessageOn Stop

```
1 msg.payload={
2   "date":"","
3   "name":"","
4   "time":""
5 }
6 return msg;
```

- Node-Red Flow to create a form




➤ Web User Interface for entering medicine details




ii) Sending the medicine name at the appropriate time:

- A function is created to compare the present time with the scheduled time.
- The function to get the present time is shown below.




```
1 var d = new Date();
2 var utc = d.getTime() + (d.getTimezoneOffset() * 60000);
3 var offset = 5.5;
4 newDate = new Date(utc + (3600000*offset));
5 var n=newDate.toISOString()
6 var date = n.slice(0,10)
7 var time = n.slice(11,16)
8 global.set('time',time)
9 msg.payload=date+" "+time
10 return msg;
11
12
```


- After obtaining the present time, cloudant node is connected to the function in order search the scheduled time by _id
- The cloudant node is connected to the IBM cloudant database.
- If the present time and scheduled time in the database matches, then the name of the medicine will be sent to the IoT device using switch node and a function following that node.



➤ Overall node flow for sending the medicine name to IoT device



iii) Complete Web UI flow:




iv) Deploying Web UI:


The debug window in Node-Red depicts the medicine details.



- The following image shows the medicine details stored in the cloudant db under the medicine database.



- The following image shows the details of the medicine in the medicine database in JSON format



- The medicine reminder operation is depicted below:

- The medicine name Dolo Cold and its intake time 17:12 and the date 17.11.2022 is stored in the database.

```
11/17/2022, 5:12:13 PM node: da03aad9e77fb7ed
msg.payload: Object
  { name: "Dolo Cold", time: "17:12",
    date: "2022-11-17" }

11/17/2022, 5:12:13 PM node: da03aad9e77fb7ed
msg.payload: Object
  { _id: "2022-11-17 17:12", name:
    "Dolo Cold" }
```

- When it's time to take the medicine, the medicine name is sent as a command to the IoT device.

```
11/17/2022, 5:12:13 PM node: da03aad9e77fb7ed
msg.payload: Object
  { _id: "2022-11-17 17:12", name:
    "Dolo Cold" }

11/17/2022, 5:12:13 PM node: c42031f5d1bb6d0cd
msg.payload: Object
  { command: "Dolo Cold" }

11/17/2022, 5:12:13 PM node: c42031f5d1bb6d0cd
msg.payload: Object
```